

Gluing micro-spheres on cantilevers in aqueous solution

Michael Glaubitz, Stephan Block, Kay-E.Gottschalk

Introduction

A new procedure was developed, which is suitable for gluing silica micro-spheres (e.g. $\varnothing = 5 \mu\text{m}$) with high accuracy on a tipless CSC12 AFM-cantilever (Mikromasch, Tallinn, Estonia). As most steps are performed in aqueous solution, the benefits are: (1) applicability even in presence of strong adhesive forces between support and micro-sphere (e.g. capillary forces) and (2) possibility to correct the micro-sphere position on the cantilever (resulting in sub- μm positional error) until the glue is cured via UV-radiation.




Fig. 1: Tipless cantilever in water with 5 μm micro-spheres

Techniques to glue micro-spheres to tipless cantilevers are commonly performed in air (for detailed information see technical report [jpk-tech-spheres-on-tip.doc](#) [1]). In this context, approaches that employ thin filaments to place the micro-sphere on the cantilever are often accompanied by a discontinuous movement / jump of the micro-sphere from the filament to the cantilever. This problem is caused by capillary forces (arising between the sphere and the cantilever) and makes it often difficult to achieve a central micro-sphere position on the cantilever beam. To increase the positional accuracy, other approaches try to catch and lift up a micro-sphere (laying on a glass-slide) directly with a glue-covered cantilever. However, these procedures become unfeasible if there are strong adhesion forces

between support and micro-sphere. Here, we present a new approach which eliminates both effects as most of the steps are performed in solution. Using this procedure it is possible to position the micro-sphere with approximately 500 nm accuracy on the cantilever beam.

Protocol

The following steps have to be performed in the sequence indicated:

- (1) At first two glass-slides have to be cleaned with the help of an RCA 1-treatment [2,3]. This treatment is recommended to avoid sphere contamination.
- (2) In the second step a thin filament or a needle is used to put a line of the UV-curable and hydrophobic adhesive NOA68 (Edmund Optics GmbH, Karlsruhe, Germany) on a clean glass-slide. It is recommended to create the line perpendicular to the orientation/long axis of the cantilever beam.
- (3) With the help of the JPK NanoWizard[®] AFM-piezos, the front of the cantilever is immersed 1 μm into the adhesive. After 1 second of indentation, the cantilever is retracted. A small drop of glue remains on the bottom side of the front of the cantilever.
- (4) 1 ml of a solution of distilled water and micro-spheres (0.006 wt. %) is pipetted on a second clean glass-slide. Because of the RCA 1-cleaning, the solution spreads over the whole slide, so that most of the micro-spheres inside the solution get separated. In the solution the micro-spheres lay on the ground, but they do not adhere. After pipetting a drop of distilled water on the cantilever, the AFM-stepper motors are used for driving the cantilever into the solution.
- (5) When catching a single micro-sphere it is important to use low tip-velocities ($\sim 1 \mu\text{m/s}$) close

to micro-sphere, because the micro-sphere could be moved away from its former position due to the induced flow.

- (6) When the glue-covered end of the cantilever beam is over the micro-sphere, the z-piezo is used to press the micro-sphere into the adhesive with a force of 5 nN for ≥ 5 s. Small adjustments of the position of the micro-sphere can be achieved via cantilever-movement with the help of the x- and y-piezo.
- (7) Now the micro-sphere weakly adheres to the tip and the cantilever can be retracted from the surface.
- (8) It is now important not to pull the cantilever out of the water, because the surface-tension could either detach the micro-sphere from the cantilever or relocate it on the cantilever. For this reason a portable UV-lamp (UVGL-25 for wavelengths of 254 and 365 nm (VWR, Darmstadt, Germany) is used for curing the adhesive in solution. Despite of the UV-absorption in water, enough radiation reaches the cantilever, so that the glue starts curing. Usually, this takes 7 minutes at 254 nm wavelength and a distance of approximately 10 cm between lamp and cantilever. After this first curing it is save to pull the cantilever out of the water. Afterwards, the curing-process has to be completed with the UV-lamp or another stronger UV-source.

central position is reached. This results reproducibly in sub- μ m positional error.

Literature

- [1] JPK technical report: Attaching microspheres to cantilevers using the NanoWizard® Life Science stage and AFM head. <http://www.jpk.com/jpk-tech-spheres-on-tip2.download.a7e6fc3981828b6eefb1cbb0557a886a.pdf>
- [2] W. Kern and J. Vossen, Eds., *Thin Film Processes*, Academic Press: New York, 1978, Ch V-1.
- [3] W. Kern, Ed., *Handbook of Semiconductor Cleaning Technology*, Noyes Publishing: Park Ridge, NJ, 1993, Ch 1

Authors adress:

Michael Glaubitz, Stephan Block, Kay-E. Gottschalk

ZIK HIKE Zentrum für Innovationskompetenz - Humorale Immunreaktionen bei kardiovaskulären Erkrankungen
Klinik und Poliklinik für HNO-Heilkunde
Fleischmannstr. 42-44
17489 Greifswald
Germany

michael.glaubitz@uni-greifswald.de




Fig. 2: Cantilever with 5 μ m micro-sphere attached at the end of the cantilever beam. Using the new deposition procedure, the location of the micro-sphere can be easily adjusted until the



HIKE

Humorale Immunreaktionen
bei kardiovaskulären Erkrankungen

ERNST MORITZ ARNDT
UNIVERSITÄT GREIFSWALD



Wissen
lockt.
Seit 1456